

For Research Use Only

PHF19 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83892-2-PBS



Basic Information

Catalog Number:

83892-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2511

GenBank Accession Number:

BC022374

GeneID (NCBI):

26147

UNIPROT ID:

Q5T6S3-2

Full Name:

PHD finger protein 19

Calculated MW:

580 aa, 63 kDa

Purification Method:

Protein A purification

CloneNo.:

241047B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00868-3, PHF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83892-6-PBS. Detection antibody: 83892-2-PBS. Standard: Ag2511 Range: 0.313-40 ng/mL