

For Research Use Only

# MPZ / PO Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83895-2-PBS



## Basic Information

Catalog Number:

83895-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0848

GenBank Accession Number:

BC006491

GeneID (NCBI):

4359

UNIPROT ID:

P25189

Full Name:

myelin protein zero

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

240948C3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP00834-1, MPZ / P0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83895-3-PBS. Detection antibody: 83895-2-PBS. Standard: Ag0848. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00834-2, MPZ / P0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83895-2-PBS. Detection antibody: 83895-1-PBS. Standard: Ag0848. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00834-3, MPZ,P0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83895-4-PBS. Detection antibody: 83895-2-PBS. Standard: Ag0848. Range: 1.563-100 ng/mL