

For Research Use Only

CD25/IL-2RA Recombinant antibody, PBS Only (Capture)

Catalog Number: 83896-2-PBS



Basic Information

Catalog Number:

83896-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC156627

GeneID (NCBI):

3559

UNIPROT ID:

P01589

Full Name:

interleukin 2 receptor, alpha

Calculated MW:

272 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

240906B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00837-2, CD25/IL-2RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83896-2-PBS. Detection antibody: 83896-1-PBS. Standard: Eg0408. Range: 2.5-40 ng/mL.