

For Research Use Only

DOT1L Recombinant antibody, PBS Only (Detector)

Catalog Number: 83904-2-PBS



Basic Information

Catalog Number:

83904-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27160

GenBank Accession Number:

NM_032482

GeneID (NCBI):

84444

UNIPROT ID:

Q8TEK3

Full Name:

DOT1-like, histone H3
methyltransferase (S. cerevisiae)

Calculated MW:

185 kDa

Purification Method:

Protein A purification

CloneNo.:

240970B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00836-1, DOT1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83904-3-PBS. Detection antibody: 83904-2-PBS. Standard: Ag27160. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83904-2-PBS against Human DOT1L was performed. The affinity constant is 2.58 nM.