

For Research Use Only

DOT1L Recombinant antibody, PBS Only (Capture)

Catalog Number: 83904-4-PBS



Basic Information

Catalog Number:

83904-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27160

GenBank Accession Number:

NM_032482

GeneID (NCBI):

84444

UNIPROT ID:

Q8TEK3

Full Name:

DOT1-like, histone H3
methyltransferase (S. cerevisiae)

Calculated MW:

185 kDa

Purification Method:

Protein A purification

CloneNo.:

240970F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00836-2, DOT1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83904-4-PBS. Detection antibody: 83904-1-PBS. Standard: Ag27160. Range: 1.56-100 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83904-4-PBS against Human DOT1L was performed. The affinity constant is 2.42 nM.