

For Research Use Only

Mouse CD83 Recombinant antibody, PBS Only (Capture)

Catalog Number:83907-1-PBS



Basic Information

Catalog Number: 83907-1-PBS	GenBank Accession Number: NM_009856.3	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 12522	CloneNo.: 241002B7
Source: Rabbit	UNIPROT ID: O88324	
Isotype: IgG	Full Name: CD83 antigen	
Immunogen Catalog Number: EG1387	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00848-2, MOUSE CD83 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83907-1-PBS. Detection antibody: 83907-4-PBS. Standard: Eg1387. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00848-3, MOUSE CD83 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83907-1-PBS. Detection antibody: 83907-3-PBS. Standard: Eg1387. Range: 0.313-40 ng/mL.