

For Research Use Only

ZFP161 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83909-3-PBS

Basic Information

Catalog Number:

83909-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16125

GenBank Accession Number:

BC110519

GeneID (NCBI):

7541

UNIPROT ID:

O43829

Full Name:

zinc finger protein 161 homolog (mouse)

Calculated MW:

449 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

240853H4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00854-2, ZFP161 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83909-3-PBS. Detection antibody: 83909-2-PBS. Standard: Ag16125. Range: 0.313-20 ng/mL