

For Research Use Only

# LOX Recombinant antibody, PBS Only (Detector)



Catalog Number: 83911-1-PBS

## Basic Information

|                                      |                                       |  |
|--------------------------------------|---------------------------------------|--|
| Catalog Number:<br>83911-1-PBS       | GenBank Accession Number:<br>BC074872 | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>4015                | CloneNo.:<br>240945D4                          |
| Source:<br>Rabbit                    | UNIPROT ID:<br>P28300                 |  |
| Isotype:<br>IgG                      | Full Name:<br>lysyl oxidase           |  |
| Immunogen Catalog Number:<br>AG12479 | Calculated MW:<br>417 aa, 47 kDa      |  |

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00857-1, LOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83911-2-PBS. Detection antibody: 83911-1-PBS. Standard: Ag12479. Range: 1.56-100 ng/mL