

For Research Use Only

CHRNA4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83915-3-PBS



Basic Information

Catalog Number:

83915-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16489

GenBank Accession Number:

BC111802

GeneID (NCBI):

1146

UNIPROT ID:

P07510

Full Name:

cholinergic receptor, nicotinic,
gamma

Calculated MW:

465 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

240986D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00853-1, CHRNG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83915-3-PBS. Detection antibody: 83915-2-PBS. Standard: Ag16489. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00853-3, CHRNG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83915-3-PBS. Detection antibody: 83915-4-PBS. Standard: Ag16489. Range: 1.25-80 ng/mL