

For Research Use Only

BUB1B Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83920-4-PBS

Basic Information

Catalog Number: 83920-4-PBS	GenBank Accession Number: BC018739	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 701	CloneNo.: 240969C9
Source: Rabbit	UNIPROT ID: O60566	
Isotype: IgG	Full Name: budding uninhibited by benzimidazoles 1 homolog beta (yeast)	
Immunogen Catalog Number: AG34182	Calculated MW: 1050 aa, 120 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00823-1, BUB1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83920-1-PBS. Detection antibody: 83920-4-PBS. Standard: Ag34182. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00823-2, BUB1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83920-4-PBS. Detection antibody: 83920-3-PBS. Standard: Ag34182. Range: 0.781-100 ng/mL.