

For Research Use Only

FOXO6 Recombinant antibody

Catalog Number: 83926-1-RR



Basic Information

Catalog Number: 83926-1-RR	GenBank Accession Number: AI361654	Purification Method: Protein A purification
Concentration: 1000 ug/ml	GeneID (NCBI): 100132074	CloneNo.: 240823D10
Source: Rabbit	Full Name: forkhead box protein O6	Recommended Dilutions: WB 1:1000-1:6000 IHC 1:200-1:800 IF-P 1:50-1:500 IF/ICC 1:50-1:500
Isotype: IgG	Calculated MW: 51 kDa Observed MW: 75 kDa	

Applications

Tested Applications: WB, IHC, IF/ICC, IF-P, ELISA	Positive Controls:
Species Specificity: human, rat	WB : HeLa cells, HepG2 cells
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	IHC : rat brain tissue,
	IF-P : rat brain tissue,
	IF/ICC : A549 cells,

Background Information

FOXO6 is a member of the murine forkhead family of transcription factors. These forkhead transcription factors may play a role in maintenance. FOXO6 mRNA is predominantly expressed in the developing brain, indicating that FOXO6 may play an important role in the nervous system.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 83926-1-RR (FOXO6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using FOXO6 antibody (83926-1-RR, Clone: 240823D10) at dilution of 1:100 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using 83926-1-RR (FOXO6 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded rat brain tissue using FOXO6 antibody (83926-1-RR, Clone: 240823D10) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).