

For Research Use Only

# CCAR1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83936-2-PBS

## Basic Information

Catalog Number:

83936-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18786

GenBank Accession Number:

BC132725

GeneID (NCBI):

55749

UNIPROT ID:

Q8IX12

Full Name:

cell division cycle and apoptosis  
regulator 1

Calculated MW:

1150 aa, 133 kDa

Purification Method:

Protein A purification

CloneNo.:

240980F6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00915-1, CCAR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83936-3-PBS. Detection antibody: 83936-2-PBS. Standard: Ag18786. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00915-2, CCAR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83936-1-PBS. Detection antibody: 83936-2-PBS. Standard: Ag18786. Range: 0.781-100 ng/mL.