

For Research Use Only

# CXXC5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83940-2-PBS

## Basic Information

Catalog Number: 83940-2-PBS	GenBank Accession Number: BC017439	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51523	CloneNo.: 24103487
Source: Rabbit	UNIPROT ID: Q7LFL8	
Isotype: IgG	Full Name: CXXC finger 5	
Immunogen Catalog Number: AG9733	Calculated MW: 322aa,33 kDa; 227aa,24 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00925-2, CXXC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83940-4-PBS. Detection antibody: 83940-2-PBS. Standard: Ag9733. Range: 0.313-40 ng/mL