

For Research Use Only

EOMES/TBR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83945-1-PBS



Basic Information

Catalog Number:

83945-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27634

GenBank Accession Number:

NM_001278182

GeneID (NCBI):

8320

UNIPROT ID:

O95936

Full Name:

eomesodermin homolog (Xenopus laevis)

Calculated MW:

73 aa

Purification Method:

Protein A purification

CloneNo.:

241036D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00931-2, EOMES/TBR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83945-4-PBS. Detection antibody: 83945-1-PBS. Standard: Ag27634. Range: 0.156-20 ng/mL.