

For Research Use Only

# EOMES/TBR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83945-4-PBS



## Basic Information

Catalog Number:

83945-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27634

GenBank Accession Number:

NM\_001278182

GeneID (NCBI):

8320

UNIPROT ID:

O95936

Full Name:

eomesodermin homolog (Xenopus  
laevis)

Calculated MW:

73 aa

Purification Method:

Protein A purification

CloneNo.:

241036G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP00931-1, EOMES/TBR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83945-4-PBS. Detection antibody: 83945-3-PBS. Standard: Ag27634. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00931-2, EOMES/TBR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83945-4-PBS. Detection antibody: 83945-1-PBS. Standard: Ag27634. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00931-3, EOMES/TBR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83945-4-PBS. Detection antibody: 83945-2-PBS. Standard: Ag27634. Range: 0.156-20 ng/mL