

For Research Use Only

HSP90 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83950-3-PBS



Basic Information

Catalog Number:

83950-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3826

GenBank Accession Number:

BC023006

GeneID (NCBI):

3320

UNIPROT ID:

P07900

Full Name:

heat shock protein 90kDa alpha
(cytosolic), class A member 1

Calculated MW:

853 aa, 90 kDa

Purification Method:

Protein A purification

CloneNo.:

241050D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00910-1, HSP90 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83950-3-PBS. Detection antibody: 83950-4-PBS. Standard: Ag3826. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00910-2, HSP90 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83950-3-PBS. Detection antibody: 83950-1-PBS. Standard: Ag3826. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83950-3-PBS against Human HSP90 was performed. The affinity constant is 1.22 nM.