

For Research Use Only

# ATG4D Recombinant antibody, PBS Only (Capture)

Catalog Number: 83951-1-PBS



## Basic Information

Catalog Number:

83951-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10177

GenBank Accession Number:

BC068992

GeneID (NCBI):

84971

UNIPROT ID:

Q86TLO

Full Name:

ATG4 autophagy related 4 homolog D  
(S. cerevisiae)

Calculated MW:

474 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

241035A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00902-2, ATG4D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83951-1-PBS. Detection antibody: 83951-3-PBS. Standard: Ag10177. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00902-3, ATG4D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83951-1-PBS. Detection antibody: 83951-4-PBS. Standard: Ag10177. Range: 0.625-40 ng/mL