

For Research Use Only

ATP6V1B1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83953-4-PBS



Basic Information

Catalog Number:

83953-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6332

GenBank Accession Number:

BC063411

GeneID (NCBI):

525

UNIPROT ID:

P15313

Full Name:

ATPase, H⁺ transporting, lysosomal
56/58kDa, V1 subunit B1

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

240955F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00904-2, ATP6V1B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83953-4-PBS. Detection antibody: 83953-1-PBS. Standard: Ag6332. Range: 0.781-100 ng/mL.