

For Research Use Only

PYROXD2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83954-2-PBS



Basic Information

Catalog Number:

83954-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14701

GenBank Accession Number:

BC006131

GeneID (NCBI):

84795

UNIPROT ID:

Q8N2H3

Full Name:

chromosome 10 open reading frame

33

Calculated MW:

581 aa, 63 kDa

Purification Method:

Protein A purification

CloneNo.:

240984F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00905-1, C10orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83954-2-PBS. Detection antibody: 83954-1-PBS. Standard: Ag14701. Range: 1.25-80 ng/mL



Cytometric bead array standard curve of MP00905-2, C10orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83954-3-PBS. Detection antibody: 83954-2-PBS. Standard: Ag14701. Range: 0.625-80 ng/mL