

For Research Use Only

# TLE3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83955-1-PBS



## Basic Information

Catalog Number:

83955-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1930

GenBank Accession Number:

BC015729

GeneID (NCBI):

7090

UNIPROT ID:

Q04726

Full Name:

transducin-like enhancer of split 3  
(E(sp1) homolog, Drosophila)

Calculated MW:

772 aa, 83 kDa

Purification Method:

Protein A purification

CloneNo.:

240680A9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00909-1, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-3-PBS. Detection antibody: 83955-1-PBS. Standard: Ag1930. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00909-2, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-2-PBS. Detection antibody: 83955-1-PBS. Standard: Ag1930. Range: 0.313-40 ng/mL