

For Research Use Only

Mouse Coagulation factor II Recombinant antibody, PBS Only

Catalog Number: 83957-4-PBS



Basic Information

Catalog Number:

83957-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0943

GenBank Accession Number:

NM_010168.3

GeneID (NCBI):

14061

UNIPROT ID:

P19221

Full Name:

coagulation factor II

Calculated MW:

70 kDa

Observed MW:

70-80 kDa

Purification Method:

Protein A purification

CloneNo.:

241013D9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

mouse, rat

Background Information

Normal blood coagulation is a complex process, involving a cascade of activation of different plasma proteins, ultimately resulting in the formation of a clot, called fibrin. Coagulation Factor II (F2), also known as prothrombin or factor II, is one of the components of this chain of plasma proteins involved in blood coagulation. Prothrombin is the precursor of thrombin, which is essential in the processes of hemostasis and thrombosis. This antibody specifically recognizes mouse Coagulation Factor II.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83957-4-RR (Coagulation factor II antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83957-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83957-4-RR against Mouse Coagulation factor II were performed. The affinity constant is 16.4 nM.