

For Research Use Only

PEA15 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83958-1-PBS



Basic Information

Catalog Number:

83958-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13782

GenBank Accession Number:

BC002426

GeneID (NCBI):

8682

UNIPROT ID:

Q15121

Full Name:

phosphoprotein enriched in astrocytes

15

Calculated MW:

130 aa, 15 kDa

Purification Method:

Protein A purification

CloneNo.:

241044A12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00916-1, PEA15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83958-2-PBS. Detection antibody: 83958-1-PBS. Standard: Ag13782. Range: 3.125-200 ng/mL.