

For Research Use Only

IL-5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83963-1-PBS



Basic Information

Catalog Number:

83963-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC066282

GeneID (NCBI):

3567

UNIPROT ID:

P05113

Full Name:

interleukin 5 (colony-stimulating factor, eosinophil)

Calculated MW:

134 aa, 15 kDa

Purification Method:

Protein A purification

CloneNo.:

240915B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00921-1, IL-5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83963-1-PBS. Detection antibody: 83963-2-PBS. Standard: Eg0380. Range: 0.313-40 ng/mL