

For Research Use Only

# IL-5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83963-2-PBS



## Basic Information

Catalog Number:

83963-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC066282

GeneID (NCBI):

3567

UNIPROT ID:

P05113

Full Name:

interleukin 5 (colony-stimulating factor, eosinophil)

Calculated MW:

134 aa, 15 kDa

Purification Method:

Protein A purification

CloneNo.:

240915F3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00921-1, IL-5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83963-1-PBS. Detection antibody: 83963-2-PBS. Standard: Eg0380. Range: 0.313-40 ng/mL