

For Research Use Only

TSC22D3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83967-1-PBS



Basic Information

Catalog Number:

83967-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3015

GenBank Accession Number:

BC018148

GeneID (NCBI):

1831

UNIPROT ID:

Q99576

Full Name:

TSC22 domain family, member 3

Calculated MW:

134 aa, 15 kDa

Purification Method:

Protein A purification

CloneNo.:

241024B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00887-2, TSC22D3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83967-1-PBS. Detection antibody: 83967-4-PBS. Standard: Ag3015. Range: 0.313-40 ng/mL