

For Research Use Only

PRMT3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83969-8-PBS



Basic Information

Catalog Number:

83969-8-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13466

GenBank Accession Number:

BC037544

GeneID (NCBI):

10196

UNIPROT ID:

O60678

Full Name:

protein arginine methyltransferase 3

Calculated MW:

531 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

241092B10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00892-4, Human PRMT3 Recombinant Matched Antibody Pair - PBS only. 83969-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag13466. 83969-4-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



A549, HepG2 and HEK-293 cell lysate were measured. The human PRMT3 concentration of detected samples was determined to be 22.38 ng/mL (based on a 1.0 mg/mL extract load) in A549 cell lysate, 9.39 ng/mL (based on a 1.2 mg/mL extract load) in HepG2 cell lysate and 45.68 ng/mL (based on a 1.2 mg/mL extract load) in HEK-293 cell lysate.