

For Research Use Only

# RAG2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83971-3-PBS



## Basic Information

Catalog Number:

83971-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2393

GenBank Accession Number:

BC022397

GeneID (NCBI):

5897

UNIPROT ID:

P55895

Full Name:

recombination activating gene 2

Calculated MW:

527 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

240855F7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00894-1, RAG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83971-4-PBS. Detection antibody: 83971-3-PBS. Standard: Ag2393. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00894-3, RAG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83971-1-PBS. Detection antibody: 83971-3-PBS. Standard: Ag2393. Range: 0.78-100 ng/mL