

For Research Use Only

RAG2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83971-5-PBS

Basic Information

Catalog Number: 83971-5-PBS	GenBank Accession Number: BC022397	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5897	CloneNo.: 240855H5
Source: Rabbit	UNIPROT ID: P55895	
Isotype: IgG	Full Name: recombination activating gene 2	
Immunogen Catalog Number: AG2393	Calculated MW: 527 aa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00894-2, RAG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83971-2-PBS. Detection antibody: 83971-5-PBS. Standard: Ag2393. Range: 0.78-100 ng/mL