

For Research Use Only

# OGG1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83973-1-PBS



## Basic Information

Catalog Number: 83973-1-PBS	GenBank Accession Number: BC000657	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 4968	CloneNo.: 240940B7
Source: Rabbit	UNIPROT ID: O15527	
Isotype: IgG	Full Name: 8-oxoguanine DNA glycosylase	
Immunogen Catalog Number: AG7204	Calculated MW: 22 kDa, 36-40 kDa, 45-57 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00871-2, OGG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83973-1-PBS. Detection antibody: 83973-3-PBS. Standard: Ag7204. Range: 0.469-60 ng/mL