

For Research Use Only

# RNGTT Recombinant antibody, PBS Only (Detector)

Catalog Number: 83974-3-PBS



## Basic Information

Catalog Number:

83974-3-PBS

Concentration:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3105

GenBank Accession Number:

BC019954

GeneID (NCBI):

8732

UNIPROT ID:

O60942

Full Name:

RNA guanylyltransferase and 5'-phosphatase

Calculated MW:

597 aa, 69 kDa

Purification Method:

Protein A purification

CloneNo.:

240927D11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Cytometric bead array standard curve of MP00899-1, RNGTT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83974-1-PBS. Detection antibody: 83974-3-PBS. Standard: Ag3105. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00899-2, RNGTT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83974-2-PBS. Detection antibody: 83974-3-PBS. Standard: Ag3105. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83974-3-PBS against Human RNGTT was performed. The affinity constant is below 1 pM.