

For Research Use Only

TIMM50 Recombinant antibody

Catalog Number: 83978-1-RR



Basic Information

Catalog Number:

83978-1-RR

Size:

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17574

GenBank Accession Number:

BC121147

GeneID (NCBI):

92609

UNIPROT ID:

Q3ZCQ8

Full Name:

translocase of inner mitochondrial
membrane 50 homolog (S. cerevisiae)

Calculated MW:

456 aa, 50 kDa

Purification Method:

Protein A purification

CloneNo.:

241097C12

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

