

For Research Use Only

PNMT Recombinant antibody, PBS Only (Capture)

Catalog Number: 83979-2-PBS



Basic Information

Catalog Number:

83979-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4098

GenBank Accession Number:

BC037246

GeneID (NCBI):

5409

UNIPROT ID:

P11086

Full Name:

phenylethanolamine N-methyltransferase

Calculated MW:

282 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

241111E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00874-1, PNMT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83979-2-PBS. Detection antibody: 83979-3-PBS. Standard: Ag4098. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00874-2, PNMT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83979-2-PBS. Detection antibody: 83979-1-PBS. Standard: Ag4098. Range: 0.625-80 ng/mL



Bilayer interferometry (BLI) kinetic assay of 83979-2-PBS against Human PNMT was performed. The affinity constant is below 1 pM.