

For Research Use Only

HSPC159 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83983-3-PBS

Basic Information

Catalog Number: 83983-3-PBS	GenBank Accession Number: BC062691	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 29094	CloneNo.: 241113G7
Source: Rabbit	UNIPROT ID: Q3ZCW2	
Isotype: IgG	Full Name: galectin-related protein	
Immunogen Catalog Number: AG10109	Calculated MW: 172 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00879-1, HSPC159 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83983-1-PBS. Detection antibody: 83983-3-PBS. Standard: Ag10109. Range: 0.78-100 ng/mL.