

For Research Use Only

Histone H1.0 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83993-2-PBS



Basic Information

Catalog Number:

83993-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9982

GenBank Accession Number:

BC000145

GeneID (NCBI):

3005

UNIPROT ID:

P07305

Full Name:

H1 histone family, member 0

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

241147D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00895-2, Histone H1.0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83993-3-PBS. Detection antibody: 83993-2-PBS. Standard: Ag9982. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00895-3, Histone H1.0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83993-2-PBS. Detection antibody: 83993-1-PBS. Standard: Ag9982. Range: 0.781-100 ng/mL