

For Research Use Only

# Mouse CTLA-4/CD152 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83995-2-PBS



## Basic Information

Catalog Number:

83995-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_009843.4

GeneID (NCBI):

12477

UNIPROT ID:

P09793

Full Name:

cytotoxic T-lymphocyte-associated protein 4

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

241100B3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00907-2, MOUSE CTLA-4/CD152 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83995-1-PBS. Detection antibody: 83995-2-PBS. Standard: Eg0632. Range: 0.156-20 ng/mL