

For Research Use Only

ARMC4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84004-1-PBS

Basic Information

Catalog Number: 84004-1-PBS	GenBank Accession Number: BC140846	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55130	CloneNo.: 241122C3
Source: Rabbit	UNIPROT ID: Q5T258	
Isotype: IgG	Full Name: armadillo repeat containing 4	
Immunogen Catalog Number: AG22340	Calculated MW: 1044 aa, 116 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00935-1, ARMC4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84004-2-PBS. Detection antibody: 84004-1-PBS. Standard: Ag22340. Range: 12.5-200 ng/mL