

For Research Use Only

RIN2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84005-2-PBS



Basic Information

Catalog Number:

84005-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35011

GenBank Accession Number:

NM_018993

GeneID (NCBI):

54453

UNIPROT ID:

Q8WYP3

Full Name:

Ras and Rab interactor 2

Calculated MW:

100 kDa

Purification Method:

Protein A purification

CloneNo.:

241233A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00936-1, RIN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84005-3-PBS. Detection antibody: 84005-2-PBS. Standard: Ag35011. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00936-3, RIN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84005-1-PBS. Detection antibody: 84005-2-PBS. Standard: Ag35011. Range: 0.313-40 ng/mL