

For Research Use Only

RIN2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84005-4-PBS

Basic Information

| | | |
|--------------------------------------|--|--|
| Catalog Number: 84005-4-PBS | GenBank Accession Number: NM_018993 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 54453 | CloneNo.: 241233C8 |
| Source: Rabbit | UNIPROT ID: Q8WYP3 | |
| Isotype: IgG | Full Name: Ras and Rab interactor 2 | |
| Immunogen Catalog Number: AG35011 | Calculated MW: 100 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00936-2, RIN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84005-4-PBS. Detection antibody: 84005-3-PBS. Standard: Ag35011. Range: 0.313-40 ng/mL