

For Research Use Only

RAB36 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84008-1-PBS



Basic Information

Catalog Number:

84008-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24921

GenBank Accession Number:

BC104991

GeneID (NCBI):

9609

UNIPROT ID:

O95755

Full Name:

RAB36, member RAS oncogene family

Calculated MW:

333 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

241083A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00933-1, RAB36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84008-1-PBS. Detection antibody: 84008-2-PBS. Standard: Ag24921. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00933-2, RAB36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84008-1-PBS. Detection antibody: 84008-3-PBS. Standard: Ag24921. Range: 0.156-20 ng/mL