

For Research Use Only

CCDC91 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84010-2-PBS

Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
84010-2-PBS	BC028682	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	55297	241159C8
Source:	UNIPROT ID:	
Rabbit	Q7Z6B0	
Isotype:	Full Name:	
IgG	coiled-coil domain containing 91	
Immunogen Catalog Number:		
AG34863		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00946-1, CCDC91 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84010-2-PBS. Detection antibody: 84010-1-PBS. Standard: Ag34863. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00946-2, CCDC91 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84010-2-PBS. Detection antibody: 84010-3-PBS. Standard: Ag34863. Range: 0.313-40 ng/mL.