

For Research Use Only

SRP54 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84014-3-PBS



Basic Information

Catalog Number:

84014-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2327

GenBank Accession Number:

BC003389

GeneID (NCBI):

6729

UNIPROT ID:

P61011

Full Name:

signal recognition particle 54kDa

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

241125G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00955-2, SRP54 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84014-2-PBS. Detection antibody: 84014-3-PBS. Standard: Ag2327. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00955-3, SRP54 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84014-1-PBS. Detection antibody: 84014-3-PBS. Standard: Ag2327. Range: 0.78-100 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84014-3-PBS against Human SRP54 was performed. The affinity constant is below 1 pM.