

For Research Use Only

IGF1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84015-2-PBS



Basic Information

Catalog Number:

84015-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29197

GenBank Accession Number:

NM_001111285.2

GeneID (NCBI):

3479

UNIPROT ID:

P05019-1

Full Name:

insulin-like growth factor 1
(somatomedin C)

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

241094G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00956-1, IGF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84015-2-PBS. Detection antibody: 84015-3-PBS. Standard: Ag29197. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00956-2, IGF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84015-1-PBS. Detection antibody: 84015-2-PBS. Standard: Ag29197. Range: 0.781-100 ng/mL.



Biolayer Interferometry (BLI) kinetic assay of 84015-2-PBS against Human IGF1 was performed. The affinity constant is 9.53 nM.