

For Research Use Only

REPS2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84016-3-PBS

Basic Information

Catalog Number:

84016-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33764

GenBank Accession Number:

NM_004726

GeneID (NCBI):

9185

UNIPROT ID:

Q8NFH8

Full Name:

RALBP1 associated Eps domain containing 2

Calculated MW:

72KD

Purification Method:

Protein A purification

CloneNo.:

241136E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00957-2, REPS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84016-3-PBS. Detection antibody: 84016-1-PBS. Standard: Ag33764. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00957-3, REPS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84016-3-PBS. Detection antibody: 84016-2-PBS. Standard: Ag33764. Range: 1.25-40 ng/mL