

For Research Use Only

TRPC3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84019-3-PBS



Basic Information

Catalog Number:

84019-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13024

GenBank Accession Number:

BC093682

GeneID (NCBI):

7222

UNIPROT ID:

Q13507

Full Name:

transient receptor potential cation
channel, subfamily C, member 3

Calculated MW:

848 aa, 97 kDa

Purification Method:

Protein A purification

CloneNo.:

241075D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00972-2, TRPC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84019-1-PBS. Detection antibody: 84019-3-PBS. Standard: Ag13024. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00972-3, TRPC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84019-2-PBS. Detection antibody: 84019-3-PBS. Standard: Ag13024. Range: 0.781-100 ng/mL.