

For Research Use Only

# TRPC3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84019-4-PBS

## Basic Information

Catalog Number:	84019-4-PBS	GenBank Accession Number:	BC093682	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	7222	CloneNo.:	241075E6
Source:	Rabbit	UNIPROT ID:	Q13507		
Isotype:	IgG	Full Name:	transient receptor potential cation channel, subfamily C, member 3		
Immunogen Catalog Number:	AG13024	Calculated MW:	848 aa, 97 kDa		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00972-1, TRPC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84019-5-PBS. Detection antibody: 84019-4-PBS. Standard: Ag13024. Range: 0.781-100 ng/mL.