

For Research Use Only

NFKB2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84022-1-PBS

Basic Information

Catalog Number:	84022-1-PBS	GenBank Accession Number:	BC002844	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	4791	CloneNo.:	241126D5
Source:	Rabbit	ENSEMBL Gene ID:	ENSG00000077150		
Isotype:	IgG	UNIPROT ID:	Q00653		
Immunogen Catalog Number:	AG0662	Full Name:	nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100)		
		Calculated MW:	97 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00950-3, NFkB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84022-1-PBS. Detection antibody: 84022-4-PBS. Standard: Ag0662. Range: Ag0662.