

For Research Use Only

ACOX3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84027-3-PBS

Basic Information

Catalog Number:

84027-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10553

GenBank Accession Number:

BC017053

GeneID (NCBI):

8310

UNIPROT ID:

O15254

Full Name:

acyl-Coenzyme A oxidase 3,
pristanoyl

Calculated MW:

624 aa, 70 kDa

Purification Method:

Protein A purification

CloneNo.:

241172G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00962-1, ACOX3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84027-3-PBS. Detection antibody: 84027-1-PBS. Standard: Ag10553. Range: 0.781-100 ng/mL.