

For Research Use Only

TMPPE Recombinant antibody, PBS Only (Capture)



Catalog Number: 84028-4-PBS

Basic Information

Catalog Number:

84028-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19556

GenBank Accession Number:

BC136873

GeneID (NCBI):

643853

UNIPROT ID:

Q6ZT21

Full Name:

transmembrane protein with metallophosphoesterase domain

Calculated MW:

453 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

241180H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00963-2, TMPPE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84028-4-PBS. Detection antibody: 84028-3-PBS. Standard: Ag19556. Range: 0.625-40 ng/mL