

For Research Use Only

ARHGAP19 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84029-1-PBS



Basic Information

Catalog Number:

84029-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34671

GenBank Accession Number:

BC114490

GeneID (NCBI):

84986

UNIPROT ID:

Q14CB8

Full Name:

Rho GTPase activating protein 19

Purification Method:

Protein A purification

CloneNo.:

241210B3

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA,
Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00937-1, Human ARHGAP19 Recombinant Matched Antibody Pair - PBS only. 84029-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag34671. 84029-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



HL-60 cell lysate was measured. The human ARHGAP19 concentration of detected samples was determined to be 11.94 ng/mL (based on a 4.3 mg/mL extract load) in HL-60 cell lysate.



Cytometric bead array standard curve of MP00937-1, ARHGAP19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84029-2-PBS. Detection antibody: 84029-1-PBS. Standard: Ag34671. Range: 1.563-200 ng/mL.