

For Research Use Only

TGM1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84032-1-PBS



Basic Information

Catalog Number:

84032-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17969

GenBank Accession Number:

BC034699

GeneID (NCBI):

7051

UNIPROT ID:

P22735

Full Name:

transglutaminase 1 (K polypeptide
epidermal type I, protein-glutamine-
gamma-glutamyltransferase)

Calculated MW:

817 aa, 90 kDa

Purification Method:

Protein A purification

CloneNo.:

241218A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00951-1, TGM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84032-2-PBS. Detection antibody: 84032-1-PBS. Standard: Ag17969. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84032-1-PBS against Human TGM1 was performed. The affinity constant is 15.9 nM.